



THE SUPERCONSTRUCT FOR CHARCOT DEFORMITY




For a Demanding Patient Population,


Use the strongest device that can be tolerated by soft tissue¹

- Beams Designed to Handle Higher Bending Forces
- Thread Design Optimized to Reduce Stress Risers
- Large Shank Beams Provide Improved Bending Strength
- Available in Sizes up to 8.5mm

AXIS STRENGTH ADVANTAGE



THREAD STRENGTH MORE MATERIAL WHERE NEEDED MOST[†]



Transcend Normal Fixation Principles

2 Maximize the Mechanical Function¹

• X-Clip acts as an "Intraosseous Anchor" to improve stability in compromised bone


Bridge Beyond the Area of Injury¹

Product Sizes		
	60-140 by 5mm lengths	
4.5mm		
	60-150 by 5mm lengths	
5.5mm*		
	70-160 by 5mm lengths	
6.5mm*	;;;;;;;;;;;;;;;;;;;;;;;;;;;;;;;;;;;;;;	
7.5mm*	TITTTTTTTTTD	
	100-170 by 10mm lengths	
8.5mm By Special Re	nuest	
	yuesi	
X-Clip (sized for each beam)	Lengths: 15, 20, 25, 30mm	
* 510(k) cleared for Charcot		




X-Clip: the X-Factor for Fusion


Overcome the Obstacle of Hardware Migration



The X-Clip is an intraosseous anchor: improves stability in compromised bone.



X-Clip for better hardware retention



X-Clip limits screw "stripout."

[†] Data on file

¹ V. James Sammarco MD, Superconstructs in the Treatment of Charcot Foot Deformity: Plantar Plating, Locked Plating and Axial Screw Fixation, *Foot Ankle Clin N AM 12 (2009) 393-407*, 399.

Indications For Use: The AXIS Charcot Fixation System in diameters of 5.5, 6.5 and 7.5mm is indicated for reconstruction procedures, non-unions and fusions of bones in the foot and ankle including the metatarsals, cuneiforms, cuboid, navicular, calcaneus and talus; specific examples include: medial and lateral column fusion resulting from neuropathic osteoarthropathy (Charcot).

Indications For Use: The 4.5/8.5 screw system is intended for fixation arthrodesis of the metatarsalcuneiform, navicular-cuneiform, metatarsal-cuboid, talonavicular, and calcaneocuboid joints.

Extremity Medical and Superconstruct are trademarks of Extremity Medical, LLC. ©Extremity Medical, LLC. 2019. All rights reserved. LBL-130-99106-EN 4/2019



300 Interpace Parkway Suite 410 Parsippany, NJ 07054 Phone: 973.588.8980 Customer Service: 888.499.0079 Fax: 888.499.0542 www.extremitymedical.com