

Omni Stable AF

Ankle Fracture System

Quick Reference Guide:
Plate and Screw Info
MKT-155-99103

Plate Offering: All Options

Plate Type	Part Number	L/R	Color	Screw Holes	Overall Length	Thickness	Proximal Screws	Compression Slot?
Lateral Fibula Plate	155-11501	Left	Blue	11 Screw Holes	68.8mm	1.5mm - 1.3mm	5 Screws	None
	155-11502	Right	Gold					
	155-11701	Left	Blue	13 Screw Holes	87.5mm		7 Screws	
	155-11702	Right	Gold					



Plate Type	Part Number	L/R	Color	Screw Holes	Overall Length	Thickness	Screws	Compression Slot?
Posterior Fibula Plate	155-12501	Left	Blue	9 Screw Holes	64.6mm	1.5mm	5 Screws	None
	155-12502	Right	Gold					
	155-12701	Left	Blue	11 Screw Holes	83.3mm		7 Screws	
	155-12702	Right	Gold					



Plate Type	Part Number	L/R	Color	Screw Holes	Overall Length	Thickness	Screws	Compression Slot?
Hook Plate	155-13301	n/a	Blue	5 Screw Holes	40.4mm	1.3mm	3 Screws	Yes
	155-13502			7 Screw Holes	59.2mm		5 Screws	



Plate Type	Part Number	L/R	Color	Screw Holes	Overall Length	Thickness	Screws	Compression Slot?
Straight Plate	155-14200	n/a	Blue	2 Screw Holes	17.5mm	1.3mm	n/a	None
	155-14400			4 Screw Holes	36.3mm			
	155-14600			6 Screw Holes	55.1mm			
	155-14800			8 Screw Holes	73.9mm			
	155-14100			10 Screw Holes	92.7mm			
	155-14120			12 Screw Holes	111.5mm			



Plate Type	Part Number	L/R	Color	Screw Holes	Overall Length	Thickness	Screws	Compression Slot?
Posterior Tibia (Y-Shape) Plate	155-21300	n/a	Blue	5 Screw Holes	36.1mm	1.5mm	3 Screws	None
	155-21500			7 Screw Holes	55.1mm		5 Screws	Yes



Omni Stable AF

Ankle Fracture System

Quick Reference Guide: Plate and Screw Info MKT-155-99103

Plate Offering (Continued)

Plate Type	Part Number	L/R	Color	Screw Holes	Overall Length	Thickness	Screws	Compression Slot?
Posterior Tibia (L-Shape) Plate	155-23301	Left	Blue	5 Screw Holes	36.1mm	1.5mm	3 Screws	None
	155-23302	Right	Gold	7 Screw Holes			5 Screws	Yes
	155-23501	Left	Blue	5 Screw Holes	55.1mm		3 Screws	None
	155-23502	Right	Gold	7 Screw Holes			5 Screws	Yes



Plate Type	Part Number	L/R	Color	Screw Holes	Overall Length	Thickness	Screws	Compression Slot?
Anterior Plate	155-24500	n/a	Blue	15 Screw Holes	65.9mm	1.5mm	5 Screws	None
	155-24700			17 Screw Holes	84.8mm		7 Screws	



Plate Type	Part Number	L/R	Color	Screw Holes	Overall Length	Thickness	Screws	Compression Slot?
Medial Tibia Plate	155-25500	n/a	Blue	9 Screw Holes	68.4mm	1.5mm - 1.3mm	5 Screws	None
	155-25700			11 Screw Holes	87.6mm		7 Screws	



Omni Stable AF Screw Offerings

Plate Screws: Solid, Fully Threaded	Major/Minor Dia	Head Type	Color	Length Offerings	Size Increments	QTY in Set	Drill Color	Screwdriver
	3.5mm / 2.7mm	Locking	Green	8mm - 60mm	2mm (< 30mm)	4 each (< 30mm) 2 each (> 30mm)	Green - 2.7mm	
		Non-Locking	Silver		5mm (> 30mm)			
	2.8mm / 2.0mm	Locking	Gold	8mm - 30mm	2mm (all)		4 each (all)	
Non-Locking		Blue						



4.0mm Screws: Solid, Fully Threaded	Major/Minor Dia	Head Type	Color	Length Offerings	Size Increments	QTY in Set	Drill Color	Screwdriver
	4.0mm / 2.9mm	Non-Locking	Dark Purple	40mm - 60mm	5mm	2 each	Purple - 2.9mm	



Cannulated Lag Screws	Major/Minor Dia	Guidewire	Color	Length Offerings	Size Increments	QTY in Set	Drill Color	Screwdriver
	3.5mm / 2.5mm	1.4mm	Magenta	20mm - 50mm	5mm	1 each	Magenta (Two Bars)	
4.0mm / 3.0mm	1.6mm	Dark Gray	20mm - 50mm	5mm	2 each	No Color Band		

